ACE XS

Q and A

www.ace-technik.com

FLOW?

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RESULTS TOO HIGH? / NO CONSISTENT VALUES?

- EN: Please wait 15 min. after drinking.
DE: Bitte warten Sie nach dem Trinken 15 Minuten.
FR: Veuillez attendre 15 minutes après avoir bu.
NL: Wacht 15 min. na het drinken.
IT: Attendere 15 minuti dopo aver bevuto.
ES: Espere 15 minutos después de beber.
-

WHY CALIBRATE?

- EN: Please calibrate the device to keep measurements precise.
DE: Bitte kalibrieren Sie das Gerät, um die Messungen präzise zu halten.
FR: Veuillez calibrer l'appareil pour que les mesures restent précises.
NL: Kalibreer het apparaat om de metingen nauwkeurig te houden.
IT: Calibrare il dispositivo per mantenere la precisione delle misurazioni.
ES: Calibre el dispositivo para que las mediciones sean precisas.

Calibration service / Kalibrierservice

Scan the QR code and download the calibration return receipt.
QR-Code scannen und Kalibrier-Rücksendeschein downloaden.

EN

DE

Happy with the product?
Leave us a review.

Breathalyzer

ACE XS

Important information about your breathalyzer – EN

Introduction

With these questions, we would like to help you better understand the topic of breathalyzers and dispel common misunderstandings. If you still have any questions, please do not hesitate to contact us. We will solve your problem quickly and easily.

Please feel free to contact us:

- **WhatsApp:** +49 1607 661 414
- **Phone:** +49 8654 588 38 0
- **Telefax:** +49 8654 / 588 38-98
- **E-mail:** info@ace-technik.com

Why does my device only display FLOW or abort the test?

The appliance is working perfectly, but it may be being used incorrectly. As soon as the device detects the air flow after approx. 0.5 seconds, it starts to beep continuously. This beeping does not indicate that you should stop, but that you should continue blowing. Continue blowing with a force similar to whistling until the device automatically ends the test with a click-clack sound.

Why do I get values that are too low?

Many people misjudge their level of alcoholization. Rules of thumb such as „2 beers, 0.5 per mille“ are misleading. How quickly alcohol is metabolized depends on many different factors, including:

- Weight
- Lung volume
- Food intake (pad)
- Gender
- Physical and mental condition
- Taking medication
- Liver performance
- Pre-existing conditions

However, the biggest factors are food intake and the time factor. A good and possibly fatty meal means that alcohol is absorbed more slowly and in smaller quantities and metabolized more quickly. In addition, the body breaks down alcohol more quickly at the beginning, when the liver is still working hard - up to 0.3 per mille per hour.

Why do I get values that are too high?

Please wait at least 15 minutes after the last sip of alcohol before carrying out the test. Otherwise, there may still be alcohol in the mouth, which will greatly distort the test results. Also, only use each mouthpiece once. Moisture residues from previous tests with higher breath alcohol concentrations can evaporate when you blow again and influence the result.

Why don't I get constant values?

Please allow 3-4 minutes between tests. A breathalyzer takes the sample after about 3-5 seconds of blowing, because the most alcohol-containing breath comes from the lower part of the lungs. The alcohol molecules settle at the bottom of the lungs and must first be blown upwards, which takes about 3-5 seconds. After a test, it takes a few minutes for the alcohol molecules to settle at the bottom of the lungs again. If you take these breaks, you will get more consistent results.

What happens if I don't have the breathalyzer calibrated regularly?

Every breathalyzer on the market is a measuring device - similar to a scale - and requires calibration/adjustment over time. Over time, the device can „drift“ slightly, which means that it loses sensitivity and tends to show lower values. If the breathalyzer is not calibrated regularly - either after 1000 tests or after 1 year - the device will become less accurate. In addition, you lose the option of extending the warranty for the sensor module from 1 to 2 years.

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ACE Handels- & Entwicklungs GmbH